

Data Sheet

50ml Centrifugal Tube

Catalog Number: 4440100N

DESCRIPTION

- Description	4440100N 50ml Centrifugal Tube
- Component Material/Resin	PP/HDPE, free from heavy metal
- Colour	Transparent ,Clear
- Dimension	See drawing below
- Sterilization	E-Beam SAL 10 ⁻⁶
- Quality Control	Each lot is sampled and inspected in accordance with Standard Operating Procedures.
- Other Information	For single use only


FEATURES

- Basic Features	DNase/RNase Free, Non-pyrogenic
- Temperature	Approx. -80 to +121°C
- Autoclavability	N/A
- Centrifugation, max, RCF	12000 x g
- Surface Treatment	N/A
- Shelf Life	3 years after date of production

PACKAGING

- Package/Material	PE
- Package Unit	25pcs/unit
- Units/Box	20 Unit/box
- Lot No.	YYYY*** (year, serial number)

DRAWING



Approved

3/04/2015