

Data Sheet

100mm Cell Culture Dish

Catalog Number: 44350300N

DESCRIPTION

- Description 4450320N 100mm Cell Culture Dish, TC Treated
- Component Material/Resin PS (Polystyrene), Free from Heavy Metal
- Colour Transparent ,Clear
- Dimension See drawing below
- Sterilization E-Beam SAL 10⁻⁶
- Quality Control Each lot is sampled and inspected in accordance with Standard Operating Procedures.
- Other Information For single use only


FEATURES

- Basic Features Non-pyrogenic
- Temperature Approx.-86 to+64°C
- Autoclavability N/A
- Centrifugation, max, RCF N/A
- Surface Treatment TC Treated
- Shelf Life 3 years after date of production

PACKAGING

- Package/Material plastic bag
- Package Unit 5 Pieces/pack
- Units/Box 60 packs/box
- Lot No. YYYY*** (year, serial number)

DRAWING



Approved 03-04-2015